



n Publication number:

**0 304 292** A2

12

## **EUROPEAN PATENT APPLICATION**

2 Application number: 88307634.1

(s) tnt. Cl.4: G 01 N 33/68

2 Date of filing: 17.08.88

30 Priority: 19.08.87 GB 8719538

(3) Date of publication of application: 22.02.89 Bulletin 89/08

Designated Contracting States:
AT BE CH DE ES FR GB GR IT LI LU NL SE

 Applicant: Farmos Group Ltd. Tengströminkatu 6 P.O. Box 425 SF-20101 Turku 10 (FI)

inventor: Risteli, Juha Suvantokatu 3 A 25 SF 90140 Oulu (FI)

> Ristell, Lella Suvantokatu 3 A 25 SF 90140 Oulu (FI)

(A) Representative: Collier, Jeremy Austin Grey et al J.A.Kemp & Co. 14, South Square Gray's Inn London WC1R 5EU (GB)

Antibody to propeptide of procollagen type III, and assay method using it. ...

(gr) A trimeric aminoterminal propeptide of type III procollagen free of proteolytic activity that could degrade the propeptide may be used to produce an antibody specific to the trimeric propeptide and which has no affinity for the monomeric degraded propeptide. This antibody may be used to assay more accurately propeptide which is a measure of rate of production of collagen, which is useful in diagnosing e.g. fibrosis.

EP 0 304 292 A2

## ANTIBODY TO PROPEPTIDE OF PROCOLLAGEN TYPE III, AND ASSAY METHOD USING IT

The present invention concerns an improved assay method for the immunological determination of the aminoterminal propeptide of type III procollagen in human serum and preparation of an antiserum suitable for use in this method.

Type III collagen is a collagen type found in several connective tissues throughout the human organism. Its proportion is characteristically high in young tissue, e.g. in the early phases of wound healing or during the initial phase of development of a fibrosis. Many diseases are associated with an increase in the amount of collagen synthesized. If the degradative mechanisms cannot compensate for the increased synthesis, accumulation of collagen in the organ takes place.

Type III collagen is synthesized as a procollagen containing propeptide extensions at both ends of the molecule. The rate of type III collegen synthesis can thus be assessed by determining the amount of propeptide liberated during the conversion of procollagen to collagen. It is known to assay the aminoterminal propeptide of type III procollagen. However, in this assay method small degradation products of the propeptide interfere, making exact measurement of the propeptide impossible and producing practical disadvantages as well as difficulties in the interpretation of the results.

It would, therefore be desirable to solve the problem of interference by the propeptide degradation products and to create a quantitative method, which is quick and simple to practise, for assaying the aminoterminal propeptide of type III procollagen in human serum.

The intact propeptide is in the form of a trimer, whereas the degradation products comprise the globular Col1 domain of the propeptide in monomeric form.

The propeptide is assayed using an antibody, generally in the form of an antiserum, raised against the propeptide. However the antiserum used has an affinity for not only the intact trimeric propeptide, but, to a lesser extent, also for the monomeric form. This is owing to the nature of the propeptide antigen used to raise the antiserum.

The propeptide antigen has been found to contain bacterial collagenase bound to its collagenous Col3 domain after in vitro liberation of the propeptide from procollagen, which is the standard method for its preparation. This Col3 domain can be digested by the enzyme in vitro at a higher temperature than the procollagen proper. Other, naturally occuring proteolytic enzymes are present in the starting material for the alternative isolation procedure (e.g. in ascitic fluid from human cancer patients). Such enzymes are capable of degrading the propeptide in vivo during the immunization process thus inducing antibodies which react with monomeric Col1 degradation products of the propeptide.

The present invention provides trimeric aminoterminal propeptide of type III procollagen free from proteolytic enzymes. Preferably the propeptide is

free from enzymes capable of degrading the propeptide to its monomeric form.

The propeptide may be purified from human ascitic fluid as is known to the biochemical experts (Biochemical Journal 1985: 232; 145-150). The isolated propeptide, which on polyacrylamide slab gel electrophoresis is homogenous and thus biochemically regarded as a single protein, may be finally subjected to low pH and is chromatographed on HPLC in acidic condition, preferably less than pH3, to release the bound enzymes applying reverse phase separation. This propeptide is then free from bacterial collagenase.

The present invention also provides an antibody raised against this propeptide, and in particular an antibody which has no affinity for the monomeric form of the propeptide.

This antibody may be raised using techniques known in the art.

The present invention also provides a method of assaying trimeric aminoterminal propeptide of type III procollagen which comprises contacting a sample to be assayed with an antibody of the invention, in the presence of a label, such that the label is bound to the propeptide-antibody complex formed, and assaying the amount of bound and/or unbound label.

The immunoassay of the invention permits the determination of the aminoterminal propeptide of type III procollagen in human serum or other body fluids without serial dilutions. This is due to the fact that the serum samples now give the same slope as the standard antigen and thus the propeptide content of an unknown sample can be directly read from the standard curve. This decreases the amount of tubes needed for the assay to one third and since there are no extra dilution steps the precision of the assay is increased and intraassay and interassay variations less than 5 % can be obtained. The new assay ensures more rapid and accurate analysis of the propeptide concentration especially in those clinical conditions (e.g. the effect of growth hormone treatment in growth hormone deficient children or in non-alcoholic liver disease) where the actual changes are not very large.

The assay may be carried out by contacting with the antibody, labelled aminoterminal propeptide of type III procollagen and the sample to be assayed, separating the propeptide antibody complex so formed from the uncomplexed material and assaying the complexed or uncomplexed label. The propeptide-antibody complex may be contacted with a second antibody which is specific to the first antibody and which may be bound to a solid support. The propeptide-antibody-antibody complex may then be separated from the uncomplexed material.

The immunological assay of the propeptide, using the reagents prepared in the way described, can be carried out e.g. with a radioactive, luminescent, fluorescent or enzymic label either in the antigen or in the antibody. Either polyclonal or monoclonal

10

20

25

30

35

50

55

60

antibodies can be used. The reference interval for adults is generally 1.7 - 4.2 micrograms/lifer serum. In children between 2 years and puberty generally 4.3 - 12.3 micrograms/liter. At puberty there is a significant increase in the serum propeptide level due to increased rate of growth.

The following examples further explain the invention

## **EXAMPLE 1**

Isolation of the aminoterminal propeptide of human type III procollagen free from proteolytic enzymes

Five to ten liters of ascitic fluid removed from cancer patients for palliative reasons are precipitated with solld (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> (40% saturation). The precipitated proteins are collected by centrifugation at 15000 x g for 30 min and dissolved in 50 mM Tris/Hc1, pH 8.6, 2 M urea containing protease inhibitors (3 milligrams/liter of phenylmethylmethylsulfonyl fluoride, N-ethylmaleimide and p-hydroxymercuribenzoate and 10 mM EDTA) and are dialyzed against this buffer. The sample is then chromatographed on a DEAE-Sephacel column (5 x 50 cm) equilibrated in this buffer. Elution is carried out with a linear gradient of NaCl (0-0.4 M NaCl, 4000 + 4000 ml). The propeptide is eluted late in the gradient, clearly after the bulk of other proteins. This material is pooled, dialyzed against 0.2 M NH4HCO3, pH 7.9, and lyophilized. The sample is dissolved in 0.2 M NH<sub>4</sub> HCO<sub>3</sub> and chromatographed on a column (1.5 x 110 cm) of Sephacryl S-300 equilibrated in this solution. The fractions containing the propeptide (near the elution position of human IgG) are pooled and dialyzed against 10 mM sodium citrate-HC1 buffer, pH 4.5. They are chromatographed on a column (1.5 x 15 cm) of DEAE-Sephacel equilibrated with this buffer. With a linear gradient of NaCl (0-0.4 M, 500 + 500 ml), the propeptide is eluted between 0.2 and 0.25 M NaCl. The propeptide obtained is homogenous on polyacrylamide slab gel electrophoresis, but still contains proteolytic enzymes attached to it. These can be effectively removed by high performance liquid chromatography applying reverse phase separation in 0.1 % trifluoroacetic acid containing 10 % 2-isopropanol and eluting the bound propeptide with increasing concentrations of 2-isopropanol (10 - 70 %, 0 - 45 min). The propeptide is eluted as a sharp peak in the first half of the gradient.

## **EXAMPLE 2**

Performance of the equilibrium type of radioimmunoassay:

Ten micrograms of the aminoterminal propeptide of type III procollagen is labelled with 1 millicurie of iodine 125 by chloramine-T (5 micrograms) and the labelled propeptide is separated from free lodine by gel filtration on a Sephacryl S-300 column (1 x 20 cm)

equilibrated in PBS-buffer. The labelled propeptide is eluted from the column as a sharp peak well before free iodine. Antiserum binding curves are prepared with 50000 radioactivity counts per minute of the labelled propeptide. The propeptide concentration in an unknown sample of serum or other body fluids is determined in the following radioimmuno inhibition assay: A pretested amount of the antiserum is incubated with the unknown sample and 50 000 counts per minute of the tracer for 2 hours at 37°. Then a solid phase second antibody against rabbit gamma globulin is added and after 15 min incubation at 20° the antigen bound in the immune complex is separated by centrifugation from the solution. The inhibition activity of the unknown sample is compared with the activity of the standard concentrations of unlabelled type III procollagen aminoterminal propeptide.

#### Claims

- 1. Trimeric aminoterminal propeptide of type III procollagen free from proteolytic enzymes.
- A propeptide as claimed in claim 1 which is free from enzymes capable of degrading the propeptide to its monomeric form.
- 3. A propertide as claimed in claim 1 or 2 free from bacterial collagenase.
- 4. A process for producing a propeptide as claimed in any one of claims 1 to 4 comprising chromatographing the propeptide under acidic conditions applying reverse phase separation.
- 5. A process according to claim 4 in which the chromatography is carried out at a pH of less than 3.
- 6. An antibody raised against a propeptide as claimed in claim 1, 2 or 3 or a propeptide produced by a process according to claim 4 or 5
- An antibody as claimed in claim 6 which has no affinity for the monomeric form of the propeptide.
- 8. A method of assaying trimeric aminoterminal propeptide of type III procollagen which comprises contacting a sample to be assayed with an antibody as claimed in claim 6 or 7, in the presence of a label such that the label is bound to the propeptide-antibody complex formed, and assaying the amount of bound and/or unbound label.
- 9. A method according to claim 8 in which labelled aminoterminal propeptide of type III procollagen and the sample to be assayed are contacted with the antibody, the propeptideantibody complex so formed is separated from the uncomplexed material and the complexed or uncomplexed label is assayed.
- 10. A method according to claim 8 or 9 in which the propeptide-antibody complex is contacted with a second antibody specific to the first antibody, and the propeptide-antibody-antibody complex is separated from the uncom-

65

5

10

15 .

plexed material.

11. A method according to claim 10 in which the second antibody is bound to a solid support.

## Claims for the following Contracting State: ES

- 1. A process for producing a trimeric aminoterminal propeptide of type III procollagen free from proteolytic enzymes said process comprising chromatographing the propeptide under acidic conditions applying reverse phase separation.
- A process according to claim 1 in which the propeptide is free from enzymes capable of degrading it to its monomeric form.
- 3. A process according to claim 1 or 2 in which the propeptide is free from bacterial collagenase.
- 4. A process according to any one of the preceding claims in which the chromatography is carried out at a pH of less than 3.
- A method of assaying trimeric aminoterminal propeptide of type III procollagen which

comprises contacting a sample to be assayed with an antibody raised against a propeptide produced by a process according to any one of the preceding claims, in the presence of a label such that the label is bound to the propeptide-antibody complex formed, and assaying the amount of bound and/or unbound label.

- 6. A method according to claim 5 in which labelled aminoterminal propeptide of type III procollagen and the sample to be assayed are contacted with the antibody, the propeptideantibody complex so formed is separated from the uncomplexed material and the complexed or uncomplexed label is assayed.
- 7. A method according to claim 5 or 6 in which the propeptide-antibody complex is contacted with a second antibody specific to the first antibody, and the propeptide-antibody-antibody complex is separated from the uncomplexed material.
- 8. A method according to claim 7 in which the second antibody is bound to a solid support.

25

20

30

35

40

45

50

55

60

Publication number:

**0 304 292** A3

# (2)

# **EUROPEAN PATENT APPLICATION**

- Application number: 88307634.1
- (51) Int. Cl.5. G01N 33/68

- ② Date of filing: 17.08.88
- (3) Priority: 19.08.87 GB 8719538
- Date of publication of application:22,02.89 Bulletin 89/08
- Designated Contracting States:

  AT BE CH DE ES FR GB GR IT LI LU NL SE
- Date of deferred publication of the search report:
   01.08.90 Bulletin 90/31
- 7) Applicant: Farmos Group Ltd.
  Tengströminkatu 6 P.O. Box 425
  SF-20101 Turku 10(FI)
- 2 Inventor: Risteli, Juha Suvantokatu 3 A 25 SF 90140 Oulu(FI) Inventor: Risteli, Leila Suvantokatu 3 A 25 SF 90140 Oulu(FI)
- Representative: Collier, Jeremy Austin Grey et al
  J.A.Kemp & Co. 14, South Square Gray's Inn
  London WC1R 5EU(GB)
- Antibody to propeptide of procollagen type III, and assay method using it.
- A trimeric aminoterminal propeptide of type III procollagen free of proteolytic activity that could degrade the propeptide may be used to produce an antibody specific to the trimeric propeptide and which has no affinity for the monomeric degraded

propeptide. This antibody may be used to assay more accurately propeptide which is a measure of rate of production of collagen, which is useful in diagnosing e.g. fibrosis.

EP 0 304 292 A3



# EUROPEAN SEARCH REPORT

Application Number

EP 88 30 7634

	DOCUMENTS CONSI	DERED TO BE RELEVAN	T	EF 66 30 70.	
Category		dication, where appropriate,	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int. Cl.4)	
D,X	BIOCHEMICAL JOURNAL vol. 232, no. 1, 1985, pages 145 - 150; O.Niem "Purification and chara N-terminal propeptide o procollagen" * page 146, column 1 - figures 1, 2 *	LONDON GB ela et al.; cterization of the f human type III	1-8	G01N33/68	
x	CLINICAL CHEMISTRY. vol. 31, no. 8, August pages 1301 - 1304; O.Ni "Radioimmunoassays for amino-terminal peptides * the whole document *	emela: type III procollagen	1-4, 6-9		
A	EP-A-4940 (MAX-PLANCK-G DER WISSENSCHAFTEN) * the whole document *	ESELLSCHAFT ZUR FÖRDERUNG	1-11		
A	EP-A-89008 (HOECHST AG)  * the whole document *		1-11	TECHNICAL FIELDS SEARCHED (Int. Cl.4)	
P.X	CLINICAL CHEMISTRY. vol. 34, no. 4, April 1 pages 715 - 718; J.Rist "Rapid equilibrium rad- amino-terminal propepti procollagen" * the whole document *	eli et al.: oimmunoassay for the	1-11	G01N C07K A61K C12P	
Ε	EP-A-289930 (HOECHST AG * the whole document *	 ) 	1-3, 6-11	 	
	The present search report has h	een drawn up for all claims			
Place of search Date of completion of the search			<del></del>	Examiner	
	BERLIN	21 MAY 1990	DE H	KOK A.J.	
X : part Y : part doc: A : tech O : non	CATEGORY OF CITED DOCUMES icularly relevant if taken alone icularly relevant if combined with and ament of the same category nological backgroundwritten disclosure rmediate document	E: earlier patent d	T: theory or principle underlying the invention E: earlier patent document, but published on, or after the filing date D: document cited in the application L: document cited for other reasons  A: member of the same patent family, corresponding document		